# beurer heart rate monitor herzfrequenz pulsuhr

german|engineering

## PM 25

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## TABLE OF CONTENTS

|   |           |
|---|-----------|
| <b>Scope of delivery</b> .....                                      | <b>2</b>  |
| <b>Important Notes</b> .....  | <b>3</b>  |
| <b>General Information for Training</b> .....                       | <b>5</b>  |
| <b>Functions of the HR monitor</b> .....                            | <b>6</b>  |
| <b>Transmission of signal and methods of Devices Measurement</b> .. | <b>7</b>  |
| <b>Getting started</b> .....  | <b>7</b>  |
| Putting on the heart rate monitor.....                              | 7         |
| Activating the heart rate monitor from the power saving mode ..     | 7         |
| Putting on the chest strap.....                                     | 8         |
| <b>General operation of the HR monitor</b> .....                    | <b>9</b>  |
| Buttons on the HR watch.....  | 9         |
| Display.....  | 10        |
| Menus.....  | 11        |
| <b>Basic settings</b> .....   | <b>13</b> |
| Overview .....  | 13        |
| Enter personal data .....   | 14        |
| Set training zone .....   | 15        |
| Set units .....   | 16        |
| <b>Watch settings</b> .....   | <b>17</b> |
| Overview .....  | 17        |
| Set time of day and date .....                                      | 17        |
| Set alarm clock.....  | 19        |
| <b>Recording of training</b> .....                                  | <b>20</b> |
| Overview .....  | 20        |
| Stop run times .....  | 21        |
| Reset recording .....   | 22        |
| <b>Results</b> .....  | <b>23</b> |
| Overview .....  | 23        |
| Display training data .....   | 23        |
| <b>Technical specifications, batteries and warranty</b> .....       | <b>25</b> |
| <b>List of abbreviations</b> .....                                  | <b>28</b> |
| <b>Index</b> .....  | <b>30</b> |

## SCOPE OF DELIVERY

|                                    |  |
|------------------------------------|--|
| Heart rate monitor with wrist band |   |
| Chest strap                        |   |
| Elastic tension strap (adjustable) |  |

In addition:

- Pouch
- Operating instructions and separate short instructions

## IMPORTANT NOTES

Please carefully read through these operating instructions, keep it for future reference and also make it available to others.

### Training

- This product is not a therapeutic device. It is a training instrument, which was developed for the measurement and representation of the human heart rate.
- With high-risk sports, please note that the use of the heart rate monitor can represent an additional source of injury.
- If in doubt or for illnesses, please consult your doctor concerning the values of your upper and lower target heart rate as well as the duration and frequency of training. In such a way, you can achieve optimal results when training.
- **WARNING:** People suffering from heart and circulatory illnesses or those with pacemakers should only use this heart rate monitor after having consulted their doctor.

### Application

- The device is only intended for private use.
- This device may only be used for the purpose it was designed for and in the manner outlined in this operating instructions. Any form of improper use can be dangerous. The manufacturer is not liable for damages that are caused by improper or incorrect use.
- The heart rate monitor is waterproof. (see chapter "Technical Specifications")
- Electromagnetic disturbances may occur. (see chapter "Technical Specifications")

## **Cleaning and care**

- From time to time, carefully clean your chest strap, elastic chest strap and heart rate monitor with a soap and water solution. Then rinse off all components with clean water. Carefully dry them with a soft towel. You can wash the elastic chest strap in the washing machine at 30° without fabric softener. The chest strap is not made for the dryer!
- Store the chest strap in a clean and dry place. Dirt impairs the elasticity and functioning of the transmitter. Sweat and moisture can keep the electrodes wet and consequently activate the transmitter, which shortens battery life.

## **Repair, Accessories and Disposal**

- Prior to use, it should be ensured that the device and accessories are free of damage. If in doubt, do not use the device and consult your sales representative or the customer service address provided.
- Repairs may only be conducted by customer service or authorized dealers. Please do not in any case attempt to repair the device!
- Use the device only in combination with the accessories provided with it.
- Avoid contact with sunscreen or the like since these can cause damage to the imprint or plastic components.
- Please dispose of the device in accordance with Waste Electric and Electronic Equipment 2002/96/EC - WEEE. For inquiries, please contact the municipal authority responsible for disposal in your area.

## GENERAL INFORMATION FOR TRAINING

This heart rate monitor serves to measure the human pulse. By means of various settings, you can support your individual training program and monitor your pulse. You should monitor your pulse during exercise because you should prevent your heart from excessive strain on the one hand and on the other in order to achieve the optimal training results. The table provides you with some information on the selection of your training zone. The maximum heart rate reduces with advancing age. Information on the training heart rate must therefore always be in relation to the maximum heart rate. The following rule of thumb is used for the identification of the maximum heart rate:

$$220 - \text{age} = \text{maximum heart rate}$$

The following example applies to a 40-year-old person:  $220 - 40 = 180$

|                              | Heart/health programme                           | Fat burning zone   | Fitness programme   | Endurance training programme                               | Anaerobic training   |
|------------------------------|--|--|---|--|--|
| <b>Max. heart rate range</b> | 50 - 60%   | 60 - 70%   | 70 - 80%  | 80 - 90%   | 90 - 100%  |
| <b>Effect</b>                | Strengthens the heart and the circulatory system | The body burns the highest percentage of calories from fat. Strengthens the heart and the circulatory system, enhances fitness | Improves the respiratory and circulatory system. Ideal to promote the basic endurance | Improves speed maintenance and increases basic speed       | Overloads specific muscles. High injury risk for athletes at amateur level, risk related to the heart in case of disease |
| <b>Suitable for</b>          | Ideal for beginners                              | Weight control and loss  | Athletes at amateur level   | Ambitious athletes at amateur level, professional athletes | Only professional athletes   |
| <b>Training</b>              | Regenerative training                            |  | Fitness training  | Endurance training   | Fitness training addressing specific zones   |

See the chapter "Basic Settings", section "Setting your training zone" on how to individually set your training zone.

## FUNCTIONS OF THE HR MONITOR

### HR functions

- Exact ECG heart rate measurement
- Transmission: analogue
- Individual training zone adjustable
- Acoustic and visual alarm for exit from the training zone
- Average heart rate (average)
- Maximum training heart rate
- Calorie consumption in Kcal (basis jogging)
- Fat burning in g/oz

### Settings

- Gender
- Weight
- Unit of weight kg/lb
- Height
- Unit of height cm/inch

### Time functions

- Time of day
- Calendar (up to 2020) and weekday
- Date
- Alarm / alarm clock
- Stop watch

## TRANSMISSION OF SIGNAL AND METHODS OF DEVICES MEASUREMENT

### Receiving range of the heart rate monitor

Your heart rate monitor picks up the heart rate signals sent by the chest strap transmitter within a range of 70 centimeters.

### Sensors on the chest strap

The entire chest strap consists of two components: the actual chest strap and an elastic tension belt. On the inside of the chest strap contacting the body, there are two rectangular, corrugated sensors. Both sensors determine your heart rate ECG exactly and transmit this to the heart rate monitor.

## GETTING STARTED

### Putting on the heart rate monitor

You can wear your heart rate monitor like a wrist watch.

### Activating the heart rate monitor from the power saving mode

You can activate your heart rate monitor by pressing any key. If the heart rate monitor has not received any signal for approx. five minutes, it switches to the power saving mode. At that point, only the current time will be displayed. Press any key in order to reactivate the watch.

## Putting on the chest strap


- Tighten the chest strap with the elastic tension belt. Adjust the length of the belt such that the belt is snug, but neither too loose nor too tight. Position the belt around the chest such that the logo is facing outwards and in the correct position directly over the breast bone. For men, the chest belt should be located directly below the pectoral muscles, for women, directly below the breast.
- Since optimal contact between the skin and the heart rate sensors can not be created immediately, it may take some time until the heart rate is measured and displayed. If necessary, alter the position of the transmitter in order to achieve optimal contact. Most of the time, sufficient contact is achieved through the accumulation of sweat beneath the chest strap. You can, however, also wet the contact points on the inside of the chest strap. For this, lift the chest strap up from the skin slightly and wet the two sensors with saliva, water or ECG gel (available at the chemistry). The contact between the skin and chest strap may not be interrupted also with greater movements of the thorax such as during deep breathing. Significant chest hair can interfere with and even prevent contact.
- Put on the chest strap a few minutes before starting in order to warm it up to body temperature and establish optimal contact.

## GENERAL OPERATION OF THE HR MONITOR

### Buttons on the HR watch

- 1  Starts and stops functions and adjusts values. Press and holding accelerates the entry.
- 2  Switches between the functions within one menu and confirms settings.
- 3  Switches between the different menus.
- 4  Pressing quickly illuminates the display for five seconds.  
Pressing and holding activates/deactivates the acoustic alarm when exiting the training zone.

## Display

- 1 Top line
- 2 Middle line
- 3 Bottom line
- 4 **Symbols**

Alarm clock is activated.

Heart rate is being received.

Training zone has been exceeded.

Training zone has been underrun.


Acoustic alarm for exit from the training zone is activated.

An explanation of the display abbreviations can be found in the list of abbreviations at the end of the operating instructions.


## Menus

You can switch menus with the  Button.

## Power saving mode


If your heart rate monitor has not received any signal for approx. five minutes, it switches to the power saving mode. In that case, only the current time will be shown. Press any key to reactivate the HR monitor.

## Main display

When changing menus, the name of the menu is first shown, then your HR monitor will automatically switch to the main display of the active menu.

From there, you have to first choose a submenu or you can directly access functions.

## Overview of the menus


| Menu     | Meaning  | Description: see chapter... |
|----------|----------|-----------------------------|
| Time     | Time     | "Watch Settings"            |
| Training | Training | "Recording of Training"     |
| Result   | Result   | "Results"                   |
| Setting  | Settings | "Basic settings"            |

## BASIC SETTINGS


### Overview

In the menu **Setting** you can:


- Enter your personal data (user). Based on this data, your HR monitor determines your calorie consumption and your fat burning during training!
- Set your training zone (limits). See chapter "General Information on Training" for explanations on how to set the training zones. Your HR monitor issues an acoustic and visual alarm during your training when the training zone is exited.
- Change the units of measurement (units) for height and weight.

Once you are in the menu **Setting**, your HR monitor automatically switches to the main display. This displays **User**.

## Enter personal data

|                   |  |   |
|-------------------|--|---|
| <p>MENU</p>       | <p>Switch to the menu <b>Setting</b>.</p>  |    |
| <p>START/STOP</p> | <p><b>Gender</b> is displayed.<br/>         START/STOP Select <b>M</b>, if you are male.<br/>         START/STOP Select <b>F</b>, if you are female.</p> |    |
| <p>OPTION/SET</p> | <p><b>Weight</b> is displayed.<br/>         START/STOP Please select your weight (20-227 kg or 44-499 lbs).</p>  |   |
| <p>OPTION/SET</p> | <p><b>Height</b> is displayed.<br/>         START/STOP Set your height (80-227 cm or 2'66"-7'6" feet).</p>   |  |
| <p>OPTION/SET</p> | <p>Return to the display <b>User</b>.</p>  |  |




## Set training zone

You can set your individual upper and lower heart rate limit as described below.


|            |   |   |
|------------|---|---|
| MENU       | Switch to the menu <b>Setting</b> .   |    |
| OPTION/SET | <b>Limits</b> is displayed.   |   |
| START/STOP | <b>Lo limit</b> is displayed.<br><input type="button" value="START/STOP"/> Set the lower limit for your training zone (40-238). |   |
| OPTION/SET | <b>Hi limit</b> is displayed.<br><input type="button" value="START/STOP"/> Set the upper limit (41-239).                        |  |
| OPTION/SET | Return to the display <b>Limits</b> .   |   |

## Set units


|                   |  |   |
|-------------------|--|---|
| <p>MENU</p>       | <p>Switch to the menu <b>Setting</b>.</p>  |    |
| <p>OPTION/SET</p> | <p><b>Limits</b> is displayed.</p>   |    |
| <p>OPTION/SET</p> | <p><b>Units</b> is displayed.</p>  |    |
| <p>START/STOP</p> | <p><b>kg/lb</b> is displayed. The weight unit blinks.<br/> <p>START/STOP Switch between the units.</p> </p>              |    |
| <p>OPTION/SET</p> | <p><b>cm/inch</b> is displayed. The unit of measurement for height blinks.<br/> <p>START/STOP Switch the units.</p> </p> |  |
| <p>OPTION/SET</p> | <p>Return to the display <b>Units</b>.</p>   |  |





# WATCH SETTINGS

## Overview

|  |   |
|--|---|
| <p>In the <b>Time</b> menu you can:</p> <ul style="list-style-type: none"> <li>■ Set the time and the date.</li> <li>■ Set the alarm clock.</li> </ul>   |  |
| <p>When you are in the <b>Time</b> menu the heart rate monitor automatically activates the main display. In this mode both the weekday and the current date are displayed. On the display: <b>Fri 12.04</b> (Friday 12th April).</p> |  |





## Set time of day and date

|                             |   |   |
|-----------------------------|---|---|
| <p>MENU</p>                 | <p>Switch to the menu <b>Time</b>.</p>  |  |
| <p>START/STOP<br/>5 sec</p> | <p><b>Set Time</b> is shown briefly. <b>Hour</b> is displayed. The hour setting blinks.<br/>START/STOP Set the hour (0-24).</p> |  |
| <p>OPTION/SET</p>           | <p><b>Minute</b> is displayed. The hour setting blinks.<br/>START/STOP Set the minute (0-59).</p>                               |  |

|                   |  |  |
|-------------------|--|--|
| <p>OPTION/SET</p> | <p><b>Year</b> is displayed. The year setting blinks.</p> <p><b>START/STOP</b> Set the year. The calendar goes up to 2020.</p> |   |
| <p>OPTION/SET</p> | <p><b>Month</b> is displayed. The monthly setting blinks.</p> <p><b>START/STOP</b> Set the month (1-12).</p>                   |   |
| <p>OPTION/SET</p> | <p><b>Dayr</b> is displayed. The day setting blinks.</p> <p><b>START/STOP</b> Set the day (1-31).</p>                          |   |
| <p>OPTION/SET</p> | <p>Back to the main display.</p>   |  |

## Set alarm clock

You can turn off the alarm clock by pressing any button. On the following day, we will be woken up at the same time.

|                             |  |   |
|-----------------------------|--|---|
| <p>MENU</p>                 | <p>Switch to the menu <b>Time</b>.</p>   |    |
| <p>OPTION/SET</p>           | <p><b>Alarm</b> is displayed.</p>  |    |
| <p>START/STOP<br/>5 sec</p> | <p><b>Set Time</b> is shown briefly.<br/><b>On Alarm</b> is displayed.<br/>START/STOP Turns the alarm clock off (Off) and on (On).</p> |   |
| <p>OPTION/SET</p>           | <p><b>Hour</b> is displayed. The hour setting blinks.<br/>START/STOP Set the hour.</p>   |  |
| <p>OPTION/SET</p>           | <p><b>Minute</b> is displayed. The hour setting blinks.<br/>START/STOP Set the minute.</p>   |  |
| <p>OPTION/SET</p>           | <p>Back to the main display.</p>   |  |

## RECORDING OF TRAINING

### Overview


In the menu **Training** you can:

- Stop your run times and thereby record your training session.
- Reset recording to zero.

The HR monitor display shows a heart rate of 138 bpm, a zeroed-out timer (00:00:00), and the word "Training" at the bottom.

Once you are in the menu **Training**, your HR monitor automatically switches to the main display. This displays the already recorded training time on the middle line (after reset: 00:00:00) and on the line below **Start**.


The HR monitor display shows a heart rate of 138 bpm, a zeroed-out timer (00:00:00), and the word "Start" at the bottom.



### Things to know about the recordings

- As soon as you have started the recording, your HR monitor starts recording your training data. You can have your training data displayed in the menu **Result**.
- Your HR monitor continually records your training data until you have reset this. This also applies when you first start training on the following day. The running time also allows you to stop your training session at any time.

## Stop run times

|                   |  |  |
|-------------------|--|--|
| <p>MENU</p>       | <p>Switch to the menu <b>Training</b>. If you have already stopped a run time, the recorded training time is displayed on the middle line.</p>             |   |
| <p>START/STOP</p> | <p>Start the recording. <b>Run</b> is displayed. Even if you change to another menu during the recording, the time continues to run in the background.</p> |   |
| <p>START/STOP</p> | <p>Stop the recording. <b>Stop</b> is displayed. You can interrupt and resume your recordings at any time by pressing the button <b>START/STOP</b>.</p>    |  |



## Reset recording

|                                    |  |  |
|------------------------------------|--|--|
|                                    | <p>You are in the menu <b>Training</b>. No training is currently being recorded. <b>Stop</b> is displayed.</p> |  <p>The LCD display shows the number '138' at the top, followed by a small icon of a person. Below that is a timer reading '00:00:03'. At the bottom, the word 'Stop' is displayed.</p>   |
| <p><b>START/STOP</b><br/>5 sec</p> | <p>Reset the recording.<br/><b>Hold Reset</b> is displayed.</p>  |  <p>The LCD display shows the words 'Hold Reset' in a large, spaced-out font.</p>   |
|                                    | <p>The display automatically switches back to the main display.</p>  |  <p>The LCD display shows the number '138' at the top, followed by a small icon of a person. Below that is a timer reading '00:00:00'. At the bottom, the word 'Start' is displayed.</p> |






# RESULTS


## Overview

|  |   |
|--|---|
| <p>In the menu <b>Result</b> you can:</p> <ul style="list-style-type: none"> <li>■ Display your training data: mean and maximum heart rate, times within and outside of the training zone as well as calorie consumption and fat burning.</li> </ul> |  |
| <p>Once you are in the menu <b>Result</b>, your HR monitor automatically switches to the main display. This displays <b>bPM AVG</b>.</p>   |  |

## Display training data

|                   |   |   |
|-------------------|---|---|
| <p>MENU</p>       | <p>Switch to the menu <b>Result</b>.</p>                          |   |
|                   | <p><b>BPM AVG</b> (average training heart rate) is displayed.</p> |  |
| <p>OPTION/SET</p> | <p><b>HRmax</b> (maximum training heart rate) is displayed.</p>   |  |

|                   |   |   |
|-------------------|---|---|
| <p>OPTION/SET</p> | <p>In [min] (time within the training zone) is displayed.</p>   |    |
| <p>OPTION/SET</p> | <p>Lo [min] (time below the training zone) is displayed.</p>  |    |
| <p>OPTION/SET</p> | <p>Hi [min] (time above the training zone) is displayed.</p>  |    |
| <p>OPTION/SET</p> | <p>[kcal] (the total calorie consumption during the recording of the training session) is displayed.</p> <p>The basis for the display of the calorie consumption and the fat burning are:</p> <ul style="list-style-type: none"> <li>■ Personal data (basic settings)</li> <li>■ the time and heart rate measured during the training session.</li> </ul> <p>The calculation is conducted based on moderately fast jogging.</p> |  |

|            |  |   |
|------------|--|---|
| OPTION/SET | <b>Fat[g]</b> (Fat burning during the entire training recording) is displayed. |  |
| OPTION/SET | Back to the main display.  |   |

## TECHNICAL SPECIFICATIONS, BATTERIES AND WARRANTY

|                     |   |
|---------------------|---|
| HR monitor          | Waterproof to 30 m (can be used in the rain)  |
| Chest strap         | Waterproof (suitable for swimming)  |
| Battery HR monitor  | 3V lithium battery, type CR2032 battery life: approx. 15 months (daily use of approx. 1 hour) |
| Battery chest strap | 3V lithium battery, type CR2032 battery life: approx. 28 months (daily use of approx. 1 hour) |


### Replacing the batteries

The battery of the HR monitor may only be replaced by a specialized dealer (any lock-maker's shop) and at the service address provided. Otherwise, the claims to warranty become void.

You can replace the battery in the chest strap on your own. The battery compartment is located on the inner side. Use a coin or the battery to be used to open the compartment. The positive pole of the battery inserted must face upwards. Carefully close the battery compartment again.

Used batteries do not belong in the household garbage. Dispose of these through your electronics dealer or your local collecting point for recoverables. You are legally required to do this. We will find this symbol on batteries containing toxic chemicals:

Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury.

## **Distortion and interference**

In the HR monitor, a receiver for signals from the chest strap transmitter is integrated. During use, it may arise that the HR monitor receives signals other than those transmitted by the chest strap.

### **How can I recognize interference?**

Interference can be recognized in that the values displayed are unrealistic. If your normal heart rate amounts to 130 bpm and it is suddenly shown to be 200, then a short-term disturbance is present. Once you have removed yourself from the vicinity of the source causing the disturbance, the values will automatically resume being accurate. Please note that distorted values can mildly influence the training result displays and overall values.

### **What kind of causes for interference are there?**

In some situations, short-range electromagnetic interference can occur. In order to avoid these, keep well away from these sources.

- Nearby power lines
- Walking through security gates in airports.
- Nearby overhead lines (tram)
- Nearby traffic lights, television sets, car motors, bike computers
- Nearby exercise equipment
- Nearby mobile phones +

## **Interference from other devices for measuring heart rates**

The HR monitor as a receiver can react to interfering stray signals such as other heart rate monitors within a vicinity of one metre. This can result in incorrect measurements particularly when training in groups. If you e.g. job next to one another, it can be beneficial to wear the HR monitor on the other wrist or to increase the distance to one's training partner.

## **Warranty**

We provide a two-year warranty from the date of purchase for material and manufacturing defects in the product.

The warranty is not valid:

- in the event of damages due to improper use
- for wearable parts such as batteries or the wrist band
- for defects that are already known to the customer at the time of purchase
- due to fault on the customer's part

The warranty does not affect the customer's legal assurances. Proof of purchase is to be provided by the customer for the assertion of a warrant claim within the period of warranty. The warranty is to be asserted within a time frame of two years from the date of purchase vis-à-vis Beurer GmbH, 89077 Ulm, Germany.

In the event of warranty, the customer is entitled to the repair of the product at our own workshops or at those authorized by us.

## LIST OF ABBREVIATIONS

Abbreviations on the display in alphabetical order:

| Abbreviation | English  |
|--------------|--|
| Alarm        | Alarm  |
| AVG          | Average heart rate during training                   |
| BPM          | Heart beats per minute                               |
| bPM AVG      | Average heart rate during training, beats per minute |
| bPM Hrmax    | Maximal heart rate during training, beats per minute |
| cm/inch      | Centimeter/Inch                                      |
| Day          | Day  |
| f            | female   |
| Fat[g]       | Fat in gram  |
| Fri          | Friday   |
| Gender       | Gender   |
| Height       | Height   |
| Hi           | Time above the training zone                         |
| Hi limit     | Highest training limit                               |
| Hour         | Hour   |
| HR           | Heart Rate   |
| HRmax        | Max. heart rate during training                      |
| In           | Time within the training zone                        |
| [Kcal]       | Calorie consumption in Kcal                          |
| kg           | Kilogram   |
| kg/lb        | Kilogram/pound                                       |
| Limits       | Training limits                                      |
| Lo           | Time below the training zone                         |
| Lo limit     | Lowest training limit                                |
| m            | male   |
| Minute       | Minute   |
| [min]        | Minute   |

| <b>Abbreviation</b> | <b>English</b>       |
|---------------------|----------------------|
| Mon                 | Monday               |
| Month               | Month                |
| Reset               | Reset values to zero |
| Result              | Result               |
| Run                 | Run time is measured |
| Sat                 | Saturday             |
| Setting             | Settings (Menu)      |
| SEt                 | Setting              |
| Start               | Start                |
| Stop                | Stop                 |
| Sun                 | Sunday               |
| Time                | Time (Menu)          |
| Thu                 | Thursday             |
| Training            | Training (Menu)      |
| Tue                 | Tuesday              |
| Units               | Units                |
| User                | User                 |
| Wed                 | Wednesday            |
| Weight              | Weight               |
| Year                | Year                 |

## INDEX

### A

Accessories 4  
Alarm Clock 10  
Application 3

### B

Basic settings 13  
Buttons 9

### C

Care 4  
Cleaning 4

### D

Display 10  
Display training data 23  
Disposal 4  
Distortion and interference 26

### G

Getting started 7

### L

List of abbreviations 28

### M

Main display 12  
Menus 11  
Methods of Device Measurement  
7

### N

Notes 3

### P

Pacemakers 3

Power saving mode 11

### R

Recording 20  
Recording of training 20  
Repair 4  
Replacing the batteries 25  
Reset 22  
Reset recording 22  
Reset stopwatch 22  
Result (Menu) 12  
Results 23

### S

Scope of delivery 2  
Set alarm clock 19  
Set date 17  
Set time of day and date 17  
Set training zone 15  
Set units 16  
Setting (Menu) 12  
Stop run times 21  
Stopwatch 21  
Symbol 10

### T

Technical specifications 25  
Time (Menu) 12  
Training 5  
Training (Menu) 12  
Transmission of Signal from the  
Devices 7

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Beurer GmbH, Söflinger Str. 218, 89077 Ulm, Germany