- |            |   |    |
|------------|---|----|
| <b>D</b>   | <b>Anhang für den behandelnden Arzt</b><br>BM 95 EKG-Stick .....      | 2  |
| <b>GB</b>  | <b>Annex for attending doctor</b><br>BM 95 ECG stick .....            | 5  |
| <b>F</b>   | <b>Annexe pour le médecin traitant</b><br>Stick ECG BM 95 .....       | 8  |
| <b>E</b>   | <b>Apéndice para el médico</b><br>Dispositivo de ECG BM 95 .....      | 11 |
| <b>I</b>   | <b>Appendice per il medico curante</b><br>BM 95 stick ECG .....       | 14 |
| <b>TR</b>  | <b>Tedaviyi yürüten doktor için ek</b><br>BM 95 EKG Çubuğu .....      | 17 |
| <b>RUS</b> | <b>Приложение для лечащего врача</b><br>Электрокардиограф BM 95 ..... | 20 |
| <b>PL</b>  | <b>Załącznik dla lekarza prowadzącego</b><br>BM 95 Element EKG .....  | 23 |

## DEUTSCH

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## ENGLISH

The ECG stick from the BM 95 range is a single-channel ECG recording device for personal use. It records the ECG in time intervals of 30 [s]. The electrodes are integrated into the device. Through their positioning, they are oriented such that the results approximately show frontal leads according to Einthoven. The following applies in this context:

Method A: approx. lead II

Method B: approx. lead III

Method C: approx. lead I

Measured electrocardiograms are recorded using the ECG stick and can then be transferred to computers on which the corresponding „Beurer CardioExpert“ software has been installed. All recordings are available to users as well as the attending doctor on screen and in printed form.


All displayed electrocardiograms are scaled by amplitude and time (recording speed of 25 [mm/s]). The ECG printout enables the user to document basic patient data, the selected lead as well as information on the current ECG recording (e.g. current complaints/medication etc.). This information can be passed on to the attending doctor to assist with the assessment of the ECG recording.

### **The described presentation of the ECG recording allows**

- (1) a doctor to subsequently provide an independent assessment of the results
- (2) for an automated, simplified assessment of the ECG results by the BM 95

The automated assessment mentioned under point (2) is not a complete or detailed clinical ECG diagnosis, but rather merely provides an initial starting point for the user due to the simplified nature of the data. The status messages from the 30 [s] – single-channel – ECG recording do not replace any medical diagnosis of any results about which you may be concerned.

## The BM 95 status messages are:

	ECG recording is OK.
	Indicates a decreased heart rate (bradycardia) of less than 55 [bpm].
	Indicates an increased heart rate (tachycardia) of greater than 100 [bpm].
	Indicates one or more pauses in the cardiac cycle, each of which is longer than 2 [s].
	<p>Indicates a disruption to the rhythm during the ECG recording.</p> <p>In general, the following disruptions to the rhythm are identified:</p> <p>Supraventricular arrhythmias (including atrial fibrillation/atrial flutter (additional display FAST)/pathological sinus arrhythmias/paroxysmal atrial tachycardias (FAST is displayed if these persist)/supraventricular extrasystoles).</p>
 	<p>Indicates a disruption to the rhythm during the ECG recording.</p> <p>Altered waveform.</p> <p>Ventricular arrhythmias (including single ventricular extrasystoles/bigeminy/trigeminy/series of ventricular extrasystoles/multifocal ventricular extrasystoles/ventricular tachycardias (additional display FAST)).</p>
	<p>Altered waveform</p> <p>Ventricular arrhythmias</p>

**The BM 95 device is not suitable for clinical emergency diagnoses.** Status messages do not inform of the presence of life-threatening disruptions to the heart rhythm such as ventricular fibrillation and/or asystole.

Status messages also do not show fundamental conduction disturbances of the heart. This applies to atrioventricular block patterns as well as complete or partial bundle branch block patterns on the left or right sides.

The primary function of the BM 95 device is to enable the recording of an electrocardiogram at home, with the aim of providing a simple means of documenting changes that could occur in the electrocardiogram in the absence of medical diagnostics. It is possible that home users may incorrectly record ECGs, e.g. due to motion artefacts, incorrect positioning of the electrodes or insufficient electrical contact (the presence of marginally smaller ECG amplitudes is signalled). As a result, any subsequent ECG diagnoses may be more difficult or even prevented.

We recommend that users undergo a medical examination if pathological status messages repeatedly occur, in particular.

**Regardless of the results delivered by the BM 95, it is always necessary for users to immediately undergo a medical examination if they experience cardiac symptoms that could indicate acute medical problems.**

Whilst the recording of BM 95 electrocardiograms can never replace a medical cardiological diagnosis, the submission of electrocardiogram results may assist with the diagnosis.

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BEURER GmbH • Söflinger Str. 218 • 89077 Ulm (Germany)