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FUNCTION KEYS ON/OFF/TARE
The ON/OFF/TARE key turns the scale on and off.
SELECT UNIT (UNIT SEL.)
Press this key to select the unit of measure you require: kg, lb or oz. These scale has a capacity of 20 kg/44.09 lb/705.4 oz, division 5 g/0.01 lb/0.1 oz.
After turning the scale on, the display shows 0.000 kg/0.00 lb/0.0 oz.
HOLD (BLOCK WEIGHT FUNCTION)
Due to the baby moving, the weight shown on the display changes all the time. While the baby is laying on the tray, press the key "HOLD".

The black point (●) blinking on the display down on the left for 3 seconds indicates that the scale is automatically calculating the average of the perceived weights. When the black point stops blinking, the display shows the weight of the baby.
The weighing remains shown on the display for some seconds; then, the scale switches off automatically.
TARE
If the baby is laid on a towel or cover, use the TARE function.
Turn the scale on with the ON/OFF/TARE key.
The display shows "8888". Wait for the display to show "0.000 kg/0.00 lb/0.0 oz" and the black point (●) down on the left, then, place the towel on the tray. Press the ON/OFF/TARE key and wait for the display to show "0.000 kg/0.00 lb/0.0 oz" and the black line (–) down on the left (about 3 seconds); then, the scale is ready for weighing.

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EN English	ELECTRONIC BABY SCALE Instructions and warranty
<p>● Never use chemical or abrasive products. ● Take care that no liquids penetrate into the scale, as it is not waterproof. Do not leave or use the scale in very damp surroundings, humidity above 85% or contact with water and other liquids could cause corrosion and effect the precision of the scale.</p>	

SYMBOLS
● "LO" indicates that the batteries are flat. Replace them.
● "0-LD" shows the weight is too much.
● The scale is ready for use when the display shows "0.000 kg/0.00 lb/0.0 oz" if the instructions have not been followed correctly, the display shows "ERR" or "-----".

TECHNICAL DATA
Range of measurement: 100 g - 20 kg
Division: 5 g
2 alkaline batteries type 1.5V AAA LR03 series
Accuracy within: +/-1% + 1 division of applied weight
Working conditions: +10°C +40°C; RH ≤85%
Storage conditions: +10°C +50°C; RH ≤85%

 i Important! Read the instructions carefully

 CE This product is meant for domestic use.
The conformity, certified by the CE-marking on the device, refers to the Directive 30/2014/EU on electromagnetic compatibility.
This product is an electronic device that was tested to assure, according to the current technical knowledge, that it does not interfere with other nearby located equipment (electromagnetic compatibility) and that it is safe if used according to the advices given in the instructions for use. In case of operation anomalies, do not keep using the device and, if necessary, refer to the manufacturer directly.
DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WeE)
The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment.
At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions.
If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m2 without having to purchase a new, similar device.
This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to its improper use of the same or of parts of the same. Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

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WARRANTY
This equipment is covered by a 2-year warranty starting from the purchase date, which has to be proven by means of the **stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.** This warranty period complies with the laws in force and applies only in case the consumer is an individual.
Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden.
The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use of the product. Use the supplied accessories only; the use of other accessories invalidate the warranty cover.
Do not open the equipment for any reason whatsoever; the opening or the tampering of the equipment, definitely invalidate the warranty. Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included. Warranty decays after 2 years from the purchase date. In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may be requested by writing at info@laica.com.

REPAIR or replacement of products covered by the warranty terms are free of charge. In case of faults, refer to the dealer; DO NOT send anything directly to LAICA.
All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product.
The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment.
Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by Laica or its dealers.

CARATERISTICHS TECHNIQUES
Gamma de mesurage: 100 g - 20 kg
Division: 5 g
2 piles alcalines type AAA de 1.5V LR03
Tolérances dans: +/-1% + 1 précision
Conditions d'utilisation: +10°C +40°C; RH ≤85%
Conditions de stockage: +10°C +50°C; RH ≤85%

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