JUMPER

Instruction Manual
TENS Therapy Device
Model: JPD-ES200

Edit date: 2019-09

Version A/0

Content

EN.....	1
DE.....	12
IT.....	24
FR.....	36
ES.....	48

Product appearance diagram

Figure-1 Appearance diagram

Main unit:


Figure-2 Front and right side for the main unit




Figure-3 Back and left side for the main unit

Accessories:

① Adhesive electrode pad

② Electrode wire

Packing list:

TENS therapy device X1

Electrode pad X 4

Electrode wire X 2

AAA 1.5V dry battery X4

Instruction manual X1

Pouch X 1

Product introduction

The TENS Therapy Device can instantaneously concentrate and spread energy of power spectrum based on a specific modulation characteristics, resulting in life-simulation of actual senses about massage, acupuncture, tapping, scraping and cupping on different body sites such as shoulders, hands, joints, abdomen, waist and legs. Frequently scheduled use of the TENS Therapy Device on meridian points can play

assistant effects to treat and alleviate various physical pains, or use of the TENS Therapy

Device for massage leads to body and mind relaxation during or after work.

Product structure and composition

It is mainly composed of TENS Therapy Device main unit, electrode connection wire, self-adhesive electrode pad and dry battery.

Intended use:

For low-frequency electric pulse treatment on human acupoints.

Scope of application:





For adjuvant treatment of pain and arthralgia such as back pain, neural paralysis and muscle pain.

Product features:

1. Five treatment modes switchable on the TENS Therapy Device;
2. Tailored pain therapy with 6 pain points (shoulder, hand, joint, abdomen, waist and leg);
3. Dual-channel treatment available for two people at the same time;
4. Adjustable intensity of treatment at any time;
5. Fall-off alarm for electrode pad, followed by automatic suspension of treatment.
6. Automatic shutdown in case of no-operation within 2 minutes except working state, more energy efficiency.
7. Exquisite and compact appearance, light weight and portability;
8. Low battery indication, alarm indication for replacing battery in time.

Install the therapy device and connect to the treatment sites:

Tear open the electrode wire transparent plastic bag and take out one electrode wire; tear open the electrode pad transparent plastic bag and take out two electrode pads; connect both heads of the electrode wire with the electrode pads, and then connect the TENS Therapy Device with the electrode wire. Remove the protective films on the electrode pads and stick the adhesive surface on the treatment sites such as the waist and other sites of the human body, and keep the electrode sheet in close contact with the skin; after use, stick the protective films back for the next use.

Operating instructions:

1. Long press on "start" button for 2s for booting;
2. Short press on "site " button to select the treatment sites (shoulders, hands, joints, abdomen, waist, leg);
3. Short press on "mode" button to select treatment modes (massage, acupuncture, tapping, scraping and cupping);

massage acupuncture knocking scraping cupping

4. Long press on "time" button for 2s to set the treatment time;
 Short press on "▲" button of A or B channel to increase the treatment time;
 long press for 2s to rapidly increase treatment time; maximum treatment time is 90min;
 Short press "▼" button of A or B channel to reduce the treatment time; long
 press for 2s to quickly reduce the treatment time. The minimum treatment time
 is 5min;
 Short press on "Start" or "time" button to exit the setting and save the current
 settings;
5. Short press on " ▲" button of A or B channel to adjust the intensity and start
 running the A channel.
 - * Note: The increment or decrement unit is 1 energy value, the minimum is 0,
 and the maximum is 60.
 - * Note: Short press on "start" button for pause of intensity output and then
 Short press on "start" button for continue running;
6. When the running time is displayed as 0, the device will automatically stop

running.

- * Note: If no button is pressed for 2 minutes, the LCD will be turned off and the system will be turned off.

7. Long press on "start" button for 2s to shutdown

Tip:

1. The treatment time is 5 minutes by default or the last set time.

Notes:

1. If you aren't feeling well during treatment, stop using it immediately.
2. During treatment, the adhesive surfaces of the electrode pads should have good contact with the skin so as to prevent sense of piercing pain due to untight contact.
3. During treatment, the slight paralysis sense may occur at treatment sites, which is the result of the coupling of the device's output current with human body, belonging to a normal circumstance.
4. Users who are treated for the first time may reduce their tolerance due to stress, so that the current intensity should not be too large. It is recommended that the intensity should be gradually increased from 0 to 60. Do not adjust the intensity too high, it is appropriate to adjust to the user's acceptable level.
5. During treatment, the intensity will be automatically reset when switching the treatment mode, and the intensity needs to be re-adjusted.

Buzzer sound specification:

Beep sound	State specification:
One short sound	Key operation
Short and continuous sound at interval of 1S	Fall-off of electrode pad

Low battery specification:

When the battery voltage is lower than $4.2 \pm 0.2V$, the battery symbol appears and flashes.

Cleaning

Cleaning of main unit: wipe the main unit and the connection cable by a soft cloth dipped with a small amount of neutral detergent. Don't infiltrate the liquid into the main unit. Don't use gasoline or volatile liquid for cleaning.

Cleaning of electrode pads: wipe the dirt on the electrode pad with a damp cloth or rinse with clear water.

Contraindications

Patients with skin diseases, malignant tumors, allergies to electrode pads, and implanted cardiac pacemaker are forbidden to use the device. Consult your doctor for other disabled situations.

Troubleshooting

Fault	Probable causes	Troubleshooting procedures
Unable to boot	Severe electricity shortage	● Replace the battery
Sense of piercing pain during	The electrode pad is not stuck tightly	● Stick electrode pad again ● Replace electrode pad
Fall-off of	<ul style="list-style-type: none">● Sweat on skin● Ineffective stickiness of	● Wipe the sweat on skin and re-attach the electrode pad
Inconsistent strength between	Disconnected main unit, electrode wires and electrode	● Make sure the connection among main unit, electrode
Power failure	<ul style="list-style-type: none">● Electricity shortage	● Replace the battery before using it
Sudden alarm pause treatment in	Fall-off of electrode pad	Re-attach the electrode pads on the treatment sites

Working conditions

Ambient temperature: 5-40°C; relative humidity: ≤80%



Storage conditions

Ambient temperature: -40°C~55°C; relative humidity: 10%~93%

Technical parameters:

1. Operating voltage: 6V DC
2. Working mode: continuous operation
3. Maximum output pulse energy: for 500 Ω of load resistance, the energy value ≤300mj
4. Pulse width: 0.150-0.200ms, error ±30%
5. Pulse frequency: 1-120Hz, error ±30%
6. Rated load resistance: 500 Ω
7. Maximum output voltage peak in open circuit: ≤500V
8. Setting range of treatment intensity: 0-60; adjustable interval: 1
9. Dimensions (L*W*H) of main unit: 121.8 mm *63.6 mm *22.2 mm
10. Electrode pad dimension: 50mm*50mm
11. Main unit weight: 86g (excluding battery)

Symbols

	Manufacturer
	Date of Manufacture
	Serials number
	Authorized representative in the European Community
	CE sign with the number of the notified body
	Waste electrical materials should be sent to a dedicated collection point for recycling.
	Type BF applied part
	Protected against solid foreign objects of 12,5 mm \varnothing and greater Protection against vertically falling water drops when ENCLOSURE tilted up to 15°
	Caution, destroy it after using
	Refer to instruction manual/booklet

ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Guidance and manufacturer's declaration - electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby
RF emissions	Class B	The device is suitable for use in all establishments other than domestic and those directly connected to the public
Harmonic emissions	Not applicable	

Voltage fluctuations/ Flicker emissions IEC 61000-3-2	Not applicable	power supply network that supplies buildings used for domestic purposes	
Guidance and manufacturer's declaration — electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an			
Immunity test	IEC 60601 T	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD)	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC	±2 kV for pow	not applicable	not applicable (For INTERNALLY POWERED ME EQUIPMENT)
Surge IEC 61000-4-5	± 1 kV line(s)	not applicable	not applicable (For INTERNALLY POWERED ME)
Voltage dips, short interruptions and voltage variations on power input lines IEC 61000-4-11	<5% (>95% dip in for 0.5 cycle 40% U _T (60% in U _T) for 5 cycles 70% U _T (30% in 25 <5% (>95% dip in U _T) for sec	not applicable	not applicable (For INTERNALLY POWERED ME EQUIPMENT)

Power r	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical
------------	------	------	--

(50Hz/60Hz) magnetic field IEC			typical commercial or hospital environment.
--------------------------------------	--	--	---

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacture's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below.

The customer or the user of device should assure that it is used in such an environment.

Immunity	IEC 61000-4-3	Compliance	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz M Hz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d = 1.167\sqrt{P}$ <p style="text-align: right;">80 MHz to 800 MHz</p> $d = 2.333\sqrt{P}$ <p style="text-align: right;">800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones
 b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V] V/m.

Recommended separation distances between portable and mobile RF communications equipment and the device.		
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.		
Maximum output power of transmitter	Separation distance according to frequency of transmitter (m)	
	$d = 1.167\sqrt{P}$ 80 MHz to 800 MHz	$d = 2.333\sqrt{P}$ 800 MHz to 2.5 GHz
.01	0.117	0.233
0.1	0.369	0.738
1	1.167	2.333
10	3.689	7.379
100	11.667	23.333

For transmitters rated at a maximum output power not listed above, There commended separation distance d in meters(m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts(W) accordable to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz. the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and refraction from structures, objects and people.

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Wellkang Ltd
Shenzhen Jumper Medical Equipment Co., Ltd.