



- ~ Both ends can be tilted, the height can be adjusted twelve high positions. When loading the stratcher up to the ambulance, the height is 62-64 cm.
- ~ The bed surface is made of waterproof foamed cushion, the height of the backrest and the knee part can be adjusted by pneumatic spring which makes the wounded feel comfortable.
- ~ The surface can be used as a carrying stretcher when it is separated from the main part.
- ~ The stretcher has the vest-style safety belts, safe and reliable.
- ~ Only one parson can operate the stretcher, sending it up to the ambulance care, the surface can be used as a seprated stretcher.
- ~ The mechanical folding legs are controlled by the handles on both sides. The stretcher must be locked on fixed device which be installed on ambulance.
- ~ It is made of high-strength aluminum alloy materials, it is light-weighted, safety-used. easy to disinfect and clean.
- ~ The size of main castor : Ø 150mm, Ø 200 mm (Inverted Cover Guardril Type)

Model	High Position (L x W x H) cm	Low Position (L x W x H) cm	Maximum Angel of the Back	Maximum Angel of the Knee Part	Self - Weight (kg)	Load Bearing (kg)	Packing		
							Quantity	Size (L x W x H) cm	G.W. (kg)
YDC-3HWF	194 x 54 x 93	194 x 54 x 34	70°	20°	42	181	1	202 x 62 x 41	49