




Beurer Connect – the perfect health monitoring system

The HealthManager represents the core system. All products labelled with the HealthManager logo come with a free version of the HealthManager PC software, app and online tool. The type of product determines the transfer method: NFC (near field communication), Bluetooth® Smart technology direct to your smartphone or the conventional method with USB to your PC. If you wish, you can save the data locally or synchronise your data with any of your systems using the web version – compatible with PC, iOS or Android.



Data accessible via web and app. Access it on the move – any time and anywhere.

The HealthManager from Beurer – the core of Beurer Connect

- **Free** We offer the innovative health monitoring system free of charge with all products indicated as HealthManager products.
- **Easy to use** Extremely easy to install and user-friendly.
- **Quick and convenient** The application is designed specifically to suit your needs – with a time-saving and convenient display.
- **At home and on the move** Your personal evaluations and data can be accessed any time and anywhere.



Good to know: key symbols



resting indicator

RESTING INDICATOR

The integral resting indicator checks whether the patient was sufficiently relaxed when the measurement was taken and thus ensures the precision of the displayed measurement result.



irregular heart rhythm

IRREGULAR HEART RHYTHM DETECTION

Warns the patient about a possible disturbance to their cardiac rhythm. Additional safety with every measurement!



beurer

THE HEALTHMANAGER


Modern health management. With the free HealthManager PC software, app and online tool, your measurements are always close to hand and you can call them up at any time.



USB
NFC

DATA TRANSFER

Transfer of the measurement values via USB and NFC technology (Near Field Communication). Transfer on an NFC-compatible smartphone at a distance of a few centimetres.)



Blood pressure monitor BM 75



german|engineering

Blood pressure monitors from Beurer – fitted with the patented resting indicator for maximum reliability of your measurement results

Upper arm blood pressure monitor BM 75



- Innovative networking between smartphone and blood pressure monitor
- NFC¹ technology (near field communication)
- Patented resting indicator² for exact measurement results
- Extra-large LCD display with white backlighting
- Average value of all stored measurements and morning and evening blood pressure for the last 7 days
- 2 x 60 memory spaces
- Classification of measurements with colour chart
- Touch-sensor buttons, blue illuminated start/stop button
- Date, time and automatic switch-off
- Optimal circumference control around arm with cuff for upper arm circumferences from 24 to 36 cm
- Incl. batteries, storage pouch, cuff holder and USB cable



It is important to check your blood pressure regularly in order to maintain and monitor early signs of high blood pressure. The **BM 75** blood pressure monitor generates the corresponding response after the measurement and guarantees maximum accuracy.

The latest generation of blood pressure measurement: The BM 75 blood pressure monitor features NFC¹ technology* (near field communication) which guarantees that measured values can be easily transferred to your smartphone.

NFC will transfer measured values to your NFC¹-compatible smartphone if you place the phone over the display. It goes without saying that the BM 75 also features the patented resting indicator² and many more extras, leaving nothing to be desired.

MEASUREMENTS CAN BE MONITORED USING BEURER HEALTHMANAGER APP, ONLINE TOOL AND PC SOFTWARE

- **Good to know** With the BM 75, you get exceptional, integrated health management in addition to the blood pressure measurement function.
- **HealthManager** You can download our online HealthManager system free of charge.
- **Ready for scan** List of compatible smartphones

The resting indicator for maximum safety and reliability

Blood pressure measurement must be taken when you are **physically** and **mentally** rested. A common error is measuring blood pressure when **the circulatory system has not been sufficiently rested**, meaning the measurement results could be an incorrect indicator of resting blood pressure.

To guarantee **accuracy** when measuring blood pressure, Beurer has developed the **world's first resting indicator**. A **symbol** informs you after the measurement whether you were **sufficiently** rested mentally and physically during the measurement.

RESTING INDICATOR DISPLAY:

Circulatory system is sufficiently **at rest during measurement** ----->

Circulatory system is **not** sufficiently **at rest** ----->

If the measurements were **not** carried out when the circulatory system was sufficiently at rest the results may differ from the "at rest" values. The measurement should be repeated after a rest and relaxation period of at least 5 minutes.

Approved for use during pregnancy*

*Clinical study: Tempestat, Institut für Medizinische Forschung, Cloppenburg, Germany



¹ NFC: Data transfer with Android version of the HealthManager app
² Patent EP 1 673 009