





## **Infrared Non-Contact Temple Thermometer**

Model: HA500 Qutie **Temple Thermometer Non Contact** START **Backlight Storage Case** 

## **Clinically Accurate**

**ASTM E1965-98:2009** 

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

- 1 second measurement
- Temple Temperature
- Object temperature
- Fever alarm
- 9 memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction























Model	Gift box	Qty per carton	Carton volume
HA500 Qutie	٧	24 pcs	0.055 cbm/ctn

Storage and Transportation Condition -25°C~55°C (-13°F~131°F) ≦95% RH, 700~1060 hPa

www.qutie-rossmax.com

OBM\_Epaper\_HA500Q\_ver1546